



Cleaning GUIDE

PAINTED STEEL

CHROME

PLASTIC

LAMINATE

STAINLESS STEEL


PROTECTIVE EXAM SCREEN



Read and understand all of the instructions in this document before cleaning your product.

PAINTED STEEL, CHROME, PLASTIC AND LAMINATE

WARNING

- 
- Always comply with facility cleaning standards and guidelines*.
 - Before using any cleaner, do a spot test on a small, inconspicuous area of your product to ensure it will not cause damage to the finish.
 - Solvents, acidic or abrasive cleaners can damage surfaces and are not recommended.
 - Wipe spill away promptly and rinse several times with water.
 - Steel wool and other abrasive cleaning pads will damage the surface.
 - Do not mix different cleaning agents to clean or disinfect.
 - Rust remover, drain cleaners and toilet bowl cleaners contain harsh chemicals and will quickly cause permanent damage.
 - Hair, textile and food dyes can cause permanent stains.

*For additional information please refer to the CDC cleaning and disinfection protocols.

RECOMMENDED CLEANING PRODUCTS

Soap and water and most common cleaners/disinfectants are excellent for most cleaning and disinfecting purposes as long as the directions on the label are followed.

The following cleaning products are acceptable when used according to directions:

- AHP Accel® TB
- Bleach Mixture (1/3 cup bleach per 1 gallon of water)
- Dispatch® Hospital Cleaner Disinfectant
- Diversey™ Oxivir® Tb Wipes
- Hydrogen Peroxide
- Isopropyl Alcohol
- Metrex™ CaviWipes™
- Oxivir® Tb
- Quaternary Cleaners
- Sani-Cloth® AF3
- Sani-Cloth® HB
- Sanimaster™ 6
- Super Sani-Cloth®


RECOMMENDED CLEANING PROCEDURES

- 1 Remove debris before application of the cleaner/disinfectant.
- 2 Apply cleaner/disinfectant according to the label. Avoid using sprays on products with electronic components.
- 3 Ensure that the solution is wiped away with a clean cloth. Do not allow the cleaner to remain on product longer than the time recommended by the label.
- 4 When solutions are NOT contained in a spray bottle make sure to appropriately dispose of the solution, cleaning cloth and cleaning supplies after each cleaning.

STAINLESS STEEL

WARNING

- Always comply with facility cleaning standards and guidelines*.
- Before using any cleaner, do a spot test on a small, inconspicuous area of your product to ensure it will not cause damage to the finish.
- DO NOT clean with detergents containing chloride.
- DO NOT USE SOLVENTS in closed spaces or while smoking.



*For additional information please refer to the CDC cleaning and disinfection protocols.

RECOMMENDED CLEANING PRODUCTS

Soap and water and most common non-abrasive cleaners/disinfectants are excellent for most cleaning and disinfecting purposes as long as the directions on the label are followed.

The following cleaning products are acceptable when used according to directions:

- Warm, clean water or soapy water
- Glass cleaner
- Non-abrasive powders; a small amount of vinegar can be added for more aggressive stains
- Organic solvents that DO NOT CONTAIN chloride (acetone, methyl alcohol and mineral spirits)
- Non-abrasive, neutral household cleaners

TIPS AND GENERAL CARE

The cleaner stainless steel can be kept while in storage, being processed or during use, the greater assurance of optimum corrosion-resistance. A neutral cleaner which is low in chloride is preferred unless the surface can be thoroughly rinsed after cleaning.

Here are a few tips on the care of stainless steel:

- Use paper or another protective wrapping on the surface of the stainless steel until processing is complete.
- Handle stainless steel with clean gloves or cloths to guard against stains or finger prints.
- Avoid the use of oily rags or greasy cloths when wiping the surface.
- Do routine cleaning of exposed surfaces.
- Where possible, after cleaning, rinse thoroughly with water.
- Even the finest cleaning powders can scratch or burnish a mill-rolled finish. On polished finishes, rubbing or wiping should be done in the direction of the polish lines, NOT across them.
- Use non-chloride or non-abrasive disinfectants, such as alcohol-based disinfectants.
- Abrasives can permanently damage some colored or highly-polished finishes. If abrasives are needed, experiment first on an inconspicuous area. Soft abrasives, such as pumice, should be used.

PROTECTIVE EXAM SCREEN

WARNING

- Always comply with facility cleaning standards and guidelines*.
- Before using any cleaner, do a spot test on a small, inconspicuous area of your product to ensure it will not cause damage to the finish.
- Protective Exam Screen, gloves, screen hole covers and screen glove fixtures must be cleaned between each patient.

*For additional information please refer to the CDC cleaning and disinfection protocols.



RECOMMENDED CLEANING PRODUCTS

The following cleaning products are acceptable when used according to directions:

- Hypochlorous Acid (HOCl)
- Bleach Mixture ($\frac{1}{3}$ cup bleach per 1 gallon of water)

RECOMMENDED CLEANING PROCEDURES

- 1 Apply cleaner/disinfectant using a spray bottle on all surfaces, and leave on for at least one minute.
- 2 Wipe away the solution with a clean, cloth or paper towel. Do not allow cleaner to remain on screen longer than the recommended time.
- 3 Replace the cloth or paper towel each time you clean.
- 4 Appropriately dispose of cleaning supplies. Multi-Use Screen Gloves may be used up to 100 times with proper cleaning.

From quality products and services to patented, industry-leading innovations, LogiQuip has been providing smart, future-focused solutions exclusively to the healthcare industry since 1992. Our knowledgeable representatives partner with you to understand your needs, then design and implement the best solution for any space in your facility. With a focus on infection control, lean processes and ROI justification, our extensive line of products and services will help make your facility more effective and efficient now and for years to come.
