

WARMING CABINETS



+ LOGIQUIP[®]
Healthcare Storage Solutions

BLANKET & FLUID WARMING CABINETS

LOGIQUIP offers a variety of Blanket and Fluid Warming Cabinets. As with all our new products, warming cabinets undergo stringent testing procedures during the manufacturing process to assure long life, the highest quality and accurate temperature settings. We offer standard sizes along with custom sizes to meet specific requirements.

Standard Features

- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2 - 6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and Visual Over-Temp Alarm
- U.L. / C.U.L. approved
- FDA Listed
- OSHPD Pre-Approved

Options

- Full view glass door
- Mobile dolly or cart
- 18" Depth
- Full recess
- Pass through chamber
- Left hand door hinge
- Keyed door lock
- Temperature recorder
- Seismic Brackets
- 220/240 V, 50/60HZ
- Temperature Data Logger
- Custom Sizes



Model # TWC24

Overall Size

30"W x 26.5"D x 78.25"T
14.00 Cubic FT Storage
1 Adjustable Shelf

Upper Chamber I.D.

26"W x 23"D x 11.25"T
6.4 Cubic FT Storage
Temperature Range 90-110 F.
Accuracy +/- 1 degree of set point

Intermediate Chamber

26"W x 23"D x 10.75"T
3.7 Cubic FT Storage
Temperature Range 90-135 F.
Accuracy +/- 1 degree of set point

Lower Chamber I.D.

26"W x 23"D x 18.50"T
3.9 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy +/- 1 degree of set point

Electrical

120V, 13.5 Amp., 50/60HZ
.65 kw/hr, 2238 btu/hr
UL/CUL Listed

(Shown on Optional Dolly)



Model # DWC24

Overall Size

30"W x 26.5"D x 74.5"T
17.2 Cubic FT Storage
3 Adjustable Shelves

Upper Chamber I.D.

26"W x 23"D x 15.25"T
5.27 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy Upper chamber
+/- 1 Degree F

Lower Chamber I.D.

26"W x 23"D x 34.5"T
11.93 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy Lower chamber
+/- 3 degrees F.

Electrical

120V, 12 Amp., 50/60HZ
.75 kw/hr, 2557 btu/hr
UL/CUL Listed
FDA Registered

(Shown on Optional Dolly)



Model # DWC2424

Overall Size

24"W x 26.5"D x 74.5"T
13.23 Cubic FT Storage
3 Adjustable Shelves

Upper Chamber I.D.

20"W x 23"D x 15.25"T
4.05 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy Upper chamber
+/- 1 Degree F.

Lower Chamber I.D.

26"W x 23"D x 34.5"T
9.18 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy Lower chamber
+/- 3 degrees F

Electrical

120V, 10 Amp., 50/60HZ
.46 kw/hr, 1570 btu/hr
UL/CUL Listed
FDA Registered

(Shown on Optional Dolly)



Model # DWC2460

Overall Size Freestanding

24"W x 20.5"D x 64.5"T
7.45 Cubic FT Storage
3 Adjustable Shelves

Upper Chamber I.D.

20"W x 17"D x 13.5"T
2.65 Cubic FT Storage
Temperature Range 90-160 F each
Chamber
Accuracy +/- 1 Degree F.

Lower Chamber I.D.

20"W x 17"D x 24.5"T
4.80 Cubic FT Storage
Temperature Range 90-160 F each
Chamber
Accuracy +/- 2 Degrees F.

Electrical

120V, 10 Amp., 50/60HZ
.49 kw/hr, 1678 btu/hr
UL/CUL Listed
FDA Registered



Model # DWC36

Overall Size

30"W x 26.5"D x 74.5"T
17.59 Cubic FT Storage
4 Adjustable Shelves

Upper Chamber I.D.

26"W x 23"D x 25"T
8.63 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy Upper Chamber
+/- 2 degrees F.

Lower Chamber I.D.

26"W x 23"D x 26"T
8.96 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy Lower Chamber
+/- 2 degrees F.

Electrical

120V, 12 Amp., 50/60HZ
.69 kw/hr, 2356 btu/hr
UL/CUL Listed
FDA Registered

(Shown on Optional Dolly)



Model # DWC60

Overall Size

30"W x 20.5"D x 64.5"T
9.73 Cubic FT Storage
3 Adjustable Shelves

Upper Chamber I.D.

26"W x 17"D x 13.5"T
3.46 Cubic FT Storage
Temperature Range 90-160 F each
Chamber
Accuracy +/- 1 Degree F.

Lower Chamber I.D.

26"W x 17"D x 24.5"T
6.27 Cubic FT Storage
Temperature Range 90-160 F each
Chamber
Accuracy +/- 2 Degrees F.

Electrical

120V, 10 Amp., 50/60HZ
.49 kw/hr, 1678 btu/hr
UL/CUL Listed
FDA Registered

WARRANTY

All warming cabinets manufactured by LOGIQUIP will carry a lifetime guarantee against product craftsmanship, 1-year labor and 2-year parts guarantee. The factory will service all units with out cost to the buyer for 1-year from shipment. After the 1-year period, replacement of a defective part (labor) will be at buyer's expense. We will exchange all defective parts at no cost to the buyer for a period of 2-years from shipment. All defective parts must be returned within 30 days to ensure proper credit.



Model # SWC1518

Overall Size
18"W x 15.5"D x 22"T

Chamber I.D.
14"W x 13"D x 13"T
1.4 Cubic FT Storage
Temperature Range 90-135 F.
Accuracy +/- 1 Degrees F.

Electrical
120V, 3.5 Amp., 50/60HZ
.39 kw/hr, 1332 btu/hr
UL/CUL Listed
FDA Registered



Model # SWC1824

Overall Size
24"W x 20.5"D x 24.5"T

Chamber I.D.
20"W x 17"D x 15.25"T
3.0 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy +/- 1 Degree F.

Electrical
120V, 3.5 Amp., 50/60HZ
.45 kw/hr, 1537 btu/hr
UL/CUL Listed
FDA Registered



Model # SWC24

Overall Size
30"W x 26.5"D x 24.5"T
1 Adjustable Shelf

Chamber I.D.
26"W x 23"D x 15.25"T
5.27 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy +/- 1 Degree F.

Electrical
120V, 7.5 Amp., 50/60HZ
.30 kw/hr, 1024 btu/hr
UL/CUL Listed
FDA Registered



Model # SWC2460

Overall Size Freestanding
24"W x 20.5"D x 64.5"T
7.45 Cubic FT Storage
3 Adjustable Shelves

Upper Chamber I.D.
20"W x 17"D x 13.5"T
2.65 Cubic FT Storage
Temperature Range 90-160 F each Chamber
Accuracy +/- 1 Degree F.

Lower Chamber I.D.
20"W x 17"D x 24.5"T
4.80 Cubic FT Storage
Temperature Range 90-160 F each Chamber
Accuracy +/- 2 Degrees F.

Electrical
120V, 7.5 Amp., 50/60HZ
.49 kw/hr, 1678 btu/hr
UL/CUL Listed
FDA Registered

(Shown on Optional Dolly)

Optional Dollies



T0007 - 26"W x 34"D x 34"T
Recommended for models – SWC24, SWC36

T0008 - 30"W x 24"D x 5"T
Recommended for all taller models (specify warmer model when ordering)

T0023 - 26"W x 20"D x 34"T
Recommended for models – SWC1518, SWC1824



Standard Keyed Temperature Lockout



Optional Door Lock



Model # SWC36

Overall Size

30"W x 26.5"D x 36"T
2 Adjustable Shelves

Chamber I.D.

26"W x 23"D x 25"T
8.63 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy +/- 2 Degrees F.

Electrical

120V, 7.5 Amp., 50/60HZ
.35 kw/hr, 1195 btu/hr
UL/CUL Listed
FDA Registered



Model # SWC60

Overall Size

30"W x 20.5"D x 64.5"T
3 Adjustable Shelves

Chamber I.D.

26"W x 17"D x 51"T
13.05 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy +/- 3 Degrees F.

Electrical

120V, 7.5 Amp., 50/60HZ
.46 kw/hr, 1570 btu/hr
UL/CUL Listed
FDA Registered

(Shown on Optional Dolly)



Model # SWC72

Overall Size

30"W x 26.5"D x 74.5"T
2 Adjustable Shelves
2 fixed Shelves

Chamber I.D.

26"W x 23"D x 61"T
21 Cubic FT Storage
Temperature Range 90-160 F.
Accuracy +/- 3 Degrees F.

Electrical

120V, 7.5Amp., 50/60HZ
.53 kw/hr, 1810 btu/hr
UL/CUL Listed
FDA Registered

(Shown on Optional Dolly)

How to order Blanket & Fluid Warming Cabinets

Options Available	18" Depth	Glass Doors	Recessed Unit	LH Door Swing, Hinges on Left	Keyed Door Lock	Temperature Recorder	220/240 V, 50/60 HZ	Pass Through Chamber
Option Extension	18	G	R	LH	DL	TR	220	P

For Example: Model # DWC24 with the 18" depth option, Glass Doors, Recessed unit, Hinges on the left, Keyed door lock, 220/240 V 50/60 HZ = Order Model # DWC24-18-G-R-LH-DL-220

If you require any further assistance please feel free to call us with any questions at 1-800-665-3760.

Standard Models

Model #	Outside Dimensions
Triple Chamber Units	
TWC24	30"W x 26.5"D x 78.25"T
TWC24-G	30"W x 26.5"D x 78.25"T
Double Chamber Units	
DWC24	30"W x 26.5"D x 74.5"T
DWC24-G	30"W x 26.5"D x 74.5"T
DWC2424	24"W x 26.5"D x 74.5"T
DWC2424-G	24"W x 26.5"D x 74.5"T
DWC2460	24"W x 20.5"D x 64.5"T
DWC2460-G	24"W x 20.5"D x 64.5"T
DWC36	30"W x 26.5"D x 74.5"T
DWC36-G	30"W x 26.5"D x 74.5"T
DWC60	30"W x 20.5"D x 64.5"T
DWC60-G	30"W x 20.5"D x 64.5"T
Single Chamber Units	
SWC1518	18"W x 15.5"D x 22"T
SWC1518-G	18"W x 15.5"D x 22"T
SWC1824	24"W x 20.5"D x 24.5"T
SWC1824-G	24"W x 20.5"D x 24.5"T
SWC24	30"W x 26.5"D x 24.5"T
SWC24-G	30"W x 26.5"D x 24.5"T
SWC2460	24"W x 20.5"D x 64.5"T
SWC2460-G	24"W x 20.5"D x 64.5"T
SWC36	30"W x 26.5"D x 36"T
SWC36-G	30"W x 26.5"D x 36"T
SWC60	30"W x 20.5"D x 64.5"T
SWC60-G	30"W x 20.5"D x 64.5"T
SWC72	30"W x 26.5"D x 74.5"T
SWC72-G	30"W x 26.5"D x 74.5"T

NOTE: Due to the vast number of sizes we have available, not all our models are shown. Please let us know your specific dimensions to get the model that will best suit your needs. **If we don't currently have the size you need, we will build a custom one for you!**



PO Box 278
 Galesburg, MI 49053
 Toll Free: 1-800-665-3760
 Fax: 1-888-965-3984
www.logiquip.com

