

# LogiCell<sup>®</sup> Scope Cabinets




Model #32636SCHG-9  
with slope top option



H.D. hinges, 5 knuckle  
hospital grade

Vent portal

**10% More Capacity than  
Competitive Cabinets!**



Drip Tray

Antibacterial lower padded wall

## LogiCell<sup>®</sup> Scope Cabinets

- Full body scope racks to protect expensive equipment bodies
- Fully padded lower cabinet back wall, antimicrobial
- Drip Tray, Removable & Cleanable
- Fully vented cabinets
- Extra capacity over competitors, 9 scopes in 36" vs 8 scopes

Model #32636SCHG-9  
with slope top option

# LogiCell® Scope Cabinets



38636SCHG9-Custom



11636SC - 9



Double Scope Depth

MODEL NO.	DIMENSIONS	DESCRIPTION
11624SC-5	18"D x 27"W x 92"H	5 Scope
11636SC-9	18"D x 36"W x 92"H	9 Scope
11636SC-18	18"D x 36"W x 92"H	14 Broncho Scopes, two level
11636SC-18DB	18"D x 36"W x 92"H	18 Full Size Scopes, 1 level dbl. deep
11642SC-11	18"D x 42"W x 92"H	11 Scope
11642SC-22DB	18"D x 42"W x 92"H	22 Full Size Scopes, 1 level dbl. deep
111824ST	18"D x 27"W	Slope Top
111836ST	18"D x 36"W	Slope Top
111842ST	18"D x 42"W	Slope Top

\* Full slope top available as accessory when specified, 96"H

## 11000 Series

- Tambour door
- Frosty white melamine
- Toe kick base, 4" recess
- Drip Tray, Removable & Cleanable
- Partial slope top \*

## 32000 Series

- Glass or solid doors
- Frosty white melamine
- Toe kick base, 4" recess
- Drip Tray, Removable & Cleanable
- Flat top \*

## 38000 Series

- Stainless cabinets
- Same sequence of sizes & options as 32000 above

MODEL NO.	DIMENSIONS	DESCRIPTION
32624SCLG-5	18"D x 27"W x 92"H	5 Scope, glass door, left hinge
32624SCRG-5	18"D x 27"W x 92"H	5 Scope, glass door, right hinge
32624SCLS-5	18"D x 27"W x 92"H	5 Scope, solid door, left hinge
32624SCRS-5	18"D x 27"W x 92"H	5 Scope, solid door, right hinge
32636SCHG-9	18"D x 36"W x 92"H	9 Scope, glass door
32636SCHS-9	18"D x 36"W x 92"H	9 Scope, Solid door
32636SCHG-18	18"D x 36"W x 92"H	18 Broncho Scopes, TwoLevel, Glass Door
32636SCHS-18	18"D x 36"W x 92"H	18 Broncho Scopes, TwoLevel, Solid Door
32636SCHG-18DB	18"D x 36"W x 92"H	18 Full Size Scopes, OneLevel, Dbl Deep, Glass Door
32636SCHS-18DB	18"D x 36"W x 92"H	18 Full Size Scopes, OneLevel, Dbl Deep, Solid Door
32642SCHG-11	18"D x 42"W x 92"H	11 Scope, Glass Door
32642SCHS-11	18"D x 42"W x 92"H	11 Scope, Solid Door
32000ST		Slope Top, Sold by the linear foot

Certified formaldehyde free option by quotation.