



CASE CART, 27"D X 46"W X 55"H- (CAT.#LQCC4655)

- stainless-steel Construction • 27" D x 46" W x 55" H • 270° Door Swing • Fully Height Adjustable Shelving - sold separately • Flat Wire, Solid, and Perforated stainless-steel Shelving Sold Separately • Four (4) 6" Casters, (2) Swivel with Brake, (1) Swivel, (1) Directional Lock • Full Wrap-Around Bumper • Suitable for High Temperature Cart Washing • Bolt-On Push Handle, Bumper Assembly, Door Hinges, Door Latches, Card Holder, and Security Hasps • Limited 5-Year Warranty

SKU: LQCC4655

Categories: [case-carts](#)

PRODUCT DESCRIPTION

LogiQuip Bolted Case Carts are easily serviceable, reducing downtime, and extending the life of your Case Cart. Designed specifically to protect sterile supplies as they move from central processing to surgical suites, LogiQuip stainless-steel Case Carts provide a safe, durable, and efficient means of transportation.

ADDITIONAL INFORMATION

Dimensions	27 × 46 × 55 in
Department	sterile-processing , surgical-services
Product Line	Case Carts
Product Category	case-carts
Product Subcategory	enclosed-carts
Material	STAINLESS-STEEL
Finish	STAINLESS-STEEL