



CASE CART OVER THE ROAD, 30" D X 43.75" W X 57" H, OVER THE ROAD CASE CART-(CAT.#OTRCC3657SS)

- All Welded, stainless-steel Construction • 30" D x 43.75" W x 57" H • 270° Door Swing • Four (4) 6" Casters, (2) Swivel with Brake, (2) Directional Lock • Door Seal • Security Breakaway Tab • Full Wrap-Around Bumper • Suitable for High Temperature Cart Washing • Limited 5 Year Warranty • Shown with Adjustable Waterfall Edge Shelves

SKU: OTRCC3657SS

Categories: [over-the-road-case-cart](#)

PRODUCT DESCRIPTION

LogiQuip's Over The Road Case Carts provide a rugged, durable and safe way to transport sterile supplies when Central Processing is off-site. All welded, stainless construction, heavy duty 6" washable casters and optional accessories including a Tow Hitch, help to create a clean and safe environment.

ADDITIONAL INFORMATION

Dimensions	30 x 43.75 x 57 in
Department	sterile-processing , surgical-services
Product Line	Case Carts
Product Category	over-the-road-case-cart
Product Subcategory	over-the-road-case-cart
Material	ALL-WELDED; STAINLESS-STEEL
Finish	STAINLESS-STEEL