



**CASE CART ROLL-OUT PERFORATED
STAINLESS STEEL SHELF, 46"-
(CAT.#LQCCS46PR)**

- Perforated stainless-steel Case Cart Shelf •
- Roll-Out Shelf • Maximum Weight Capacity: 200 lbs. • Compatible with LQCC4640, LQCC4655

SKU: LQCCS46PR

Categories: [case-carts-components](#)

PRODUCT DESCRIPTION

Roll-Out Perforated stainless-steel Shelves are designed to be durable. These shelves can be placed anywhere along the track on the inside of the case cart to maximize capacity. With smooth glide rails, you will have easy access to stored items.

ADDITIONAL INFORMATION

Dimensions	23 × 35.75 × 1 in
Department	sterile-processing , surgical-services
Product Line	Case Carts
Product Category	case-carts-components
Product Subcategory	case-cart-shelves
Material	STAINLESS-STEEL
Finish	STAINLESS-STEEL

