



CASE CART ROLL-OUT SOLID STAINLESS STEEL SHELF, 36"-(CAT.#LQCCS36SR)

- Solid stainless-steel Case Cart Shelf • Roll-Out Shelf • Maximum Weight Capacity: 200 lbs. • Compatible with LQCC3630, LQCC3640, and LQCC3655


SKU: LQCCS36SR

Categories: [case-carts-components](#)

PRODUCT DESCRIPTION

Roll-Out Solid stainless-steel Shelves are designed to be durable. These shelves can be placed anywhere along the track on the inside of the case cart to maximize capacity. With smooth glide rails, you will have easy access to stored items.

ADDITIONAL INFORMATION



Dimensions	23 × 26 × 1 in
Department	sterile-processing , surgical-services
Product Line	Case Carts
Product Category	case-carts-components
Product Subcategory	case-cart-shelves
Material	STAINLESS-STEEL
Finish	STAINLESS-STEEL