



**CASE CART ROLL-OUT WIRE SHELF, 36"-
(CAT.#LQCCS36WR)**

- Wire Case Cart Shelf • Roll-Out Shelf •
- Maximum Weight Capacity: 150 lbs. •
- Compatible with LQCC3630, LQCC3640, and LQCC3655

SKU: LQCCS36WR

Categories: [case-carts-components](#)

PRODUCT DESCRIPTION

Roll-Out Wire Shelves are designed to be durable. This rounded edge helps to ensure sterile packs are not damaged when moved in and out of the case cart.

ADDITIONAL INFORMATION

Dimensions	23 × 26 × 1 in
Department	sterile-processing , surgical-services
Product Line	Case Carts
Product Category	case-carts-components
Product Subcategory	case-cart-shelves
Material	STAINLESS-STEEL
Finish	STAINLESS-STEEL

