



CASE CART TOP EXTENSION SHELF- (CAT.#CC-TSE)

- Compatible our Bolted Case Cart Line •
- stainless-steel Construction

SKU: CC-TSE

Categories: [case-carts-accessories](#)

PRODUCT DESCRIPTION

This stainless-steel shelf can be added to the side of your Case Cart for added work space. Simply flip-up and lock in place.

ADDITIONAL INFORMATION

Dimensions	24 × 16 in
Department	sterile-processing , surgical-services
Product Line	Case Carts
Product Category	case-carts-accessories
Product Subcategory	case-carts-accessories
Material	STAINLESS-STEEL
Finish	STAINLESS-STEEL

