



Product Image  
Coming Soon

**[-\(CAT.#42490SCHG-10PTD-TSRF\)](#)**

- HEPA Power Vent with Channel Purge •
- ANSI/AAMI ST91: 2021, SGNA, and AORN
- Compliant • RFID Reader • 304 Stainless Steel
- Construction • 10 Durable, Easy-to-Clean Large
- Diameter Scope Holders • LED Lighting •
- Removable Drip Tray • Maintenance Plans
- Available • Limited 5-Year Warranty

---

**SKU:** 42490SCHG-10PTD-TSRF

**Categories:** [dri-stor-scope-cabinets](#)

## PRODUCT DESCRIPTION

HEPA Power Vent with Channel Purge Technology comply with ANSI/AAMI ST91:2021 Guidelines, emphasizing the importance of drying the internal channels of flexible and semi-rigid endoscopes by forcing HEPA Filtered Air through the tubing channels assigned to individual Scopes for 30-minute cycles. Airflow is discontinued for 20 minutes before a new cycle begins The Filtered Air quality is to 0.3 micron at 99.97% efficiency. In addition, this unit has a RFID reader for controlled access.

## ADDITIONAL INFORMATION

<b>Department</b>	<a href="#">cath-lab-ir-endoscopy</a>
<b>Product Line</b>	<a href="#">Scope Cabinets</a>
<b>Product Category</b>	<a href="#">dri-stor-scope-cabinets</a>
<b>Product Subcategory</b>	<a href="#">dri-stor-10-scopes</a>
<b>Material</b>	<a href="#">STAINLESS-STEEL</a>