



Product Image
Coming Soon

[-\(CAT.#43690SCHG-20PTD-TSRF\)](#)

- HEPA Power Vent with Channel Purge •
- ANSI/AAMI ST91: 2021, SGNA, and AORN
- Compliant • RFID Reader • 304 Stainless Steel
- Construction • 20 Durable, Easy-to-Clean Large
- Diameter Scope Holders • LED Lighting •
- Removable Drip Tray • Maintenance Plans
- Available • Limited 5-Year Warranty

SKU: 43690SCHG-20PTD-TSRF

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION

HEPA Power Vent with Channel Purge Technology comply with ANSI/AAMI ST91:2021 Guidelines, emphasizing the importance of drying the internal channels of flexible and semi-rigid endoscopes by forcing HEPA Filtered Air through the tubing channels assigned to individual Scopes for 30-minute cycles. Airflow is discontinued for 20 minutes before a new cycle begins The Filtered Air quality is to 0.3 micron at 99.97% efficiency. In addition, this unit has a RFID reader for controlled access.

ADDITIONAL INFORMATION

| | |
|----------------------------|---|
| Department | cath-lab-ir-endoscopy |
| Product Line | Scope Cabinets |
| Product Category | dri-stor-scope-cabinets |
| Product Subcategory | dri-stor-20-scopes |
| Material | STAINLESS-STEEL |