



[-\(CAT.#43890SCHG-20PSD-ADD-TS\)](#)

- HEPA Power Vent • ANSI/AAMI ST91: 2021, SGNA, and AORN Compliant • Equipped with EndoTrack • 304 Stainless Steel Construction • 20 Durable, Easy-to-Clean Large Diameter Scope Holders • LED Lighting • Removable Drip Tray • Maintenance Plans Available • Limited 5-Year Warranty

SKU: 43890SCHG-20PSD-ADD-TS

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION

HEPA Power Vent with Channel Purge Technology comply with ANSI/AAMI ST91:2021 Guidelines, emphasizing the importance of drying the internal channels of flexible and semi-rigid endoscopes by forcing HEPA Filtered Air through the tubing channels assigned to individual Scopes for 30-minute cycles. Airflow is Dri-Stor Top Mount. Scope Cabinets are equipped with a HEPA Power Vent that circulates HEPA Filtered Air throughout the Cabinet for infection-control conscious drying and storage. The Filtered Air quality is to 0.3 micron at 99.97% efficiency. In addition, this unit is equipped with EndoTrack. EndoTrack allows you to easily manage and track scopes and probes through an intuitive interface. Utilizing barcode or RFID technology, EndoTrack provides controlled access and security to your cabinet as well as the ability to track hang time.

ADDITIONAL INFORMATION

Department	cath-lab-ir-endoscopy
Product Line	Scope Cabinets
Product Category	dri-stor-scope-cabinets
Product Subcategory	dri-stor-20-scopes
Material	STAINLESS-STEEL