



### **-(CAT.#43890SCHG-20PSD-ADD)**

- Compatible with LogiQuip Dri-Stor HEPA Channel Purge Side Mount Cabinets • HEPA Power Vent with Channel Purge • ANSI/AAMI ST91: 2021, SGNA, and AORN Compliant • 304 Stainless Steel Construction • 20 Durable, Easy-to-Clean Large Diameter Scope Holders • LED Lighting • Removable Drip Tray • Maintenance Plans Available • Limited 5-Year Warranty

---

**SKU:** 43890SCHG-20PSD-ADD

**Categories:** [dri-stor-scope-cabinets](#)

## **PRODUCT DESCRIPTION**

Expand the capacity of your Dri-Stor Side Mount Scope Cabinet with this Add-On unit. Compatible with our Dri-Stor Side Mount Cabinets. HEPA Power Vent with Channel Purge Technology comply with ANSI/AAMI ST91:2021 Guidelines, emphasizing the importance of drying the internal channels of flexible and semi-rigid endoscopes by forcing HEPA Filtered Air through the tubing channels assigned to individual Scopes for 30-minute cycles. Airflow is discontinued for 20 minutes before a new cycle begins. The Filtered Air quality is to 0.3 micron at 99.97% efficiency.

## **ADDITIONAL INFORMATION**

<b>Department</b>	<a href="#">cath-lab-ir-endoscopy</a>
<b>Product Line</b>	<a href="#">Scope Cabinets</a>
<b>Product Category</b>	<a href="#">dri-stor-scope-cabinets</a>
<b>Product Subcategory</b>	<a href="#">dri-stor-20-scopes</a>
<b>Material</b>	<a href="#">STAINLESS-STEEL</a>