



DECONTAMINATION SINK, ADJUSTABLE HEIGHT, 120" WIDE, TRIPLE BASIN- (CAT.#LQDSA3012035)

- 35" - 42" Adjustable Height, Electronic Control
- 120" Wide Triple Basin Sink with Etched Gallon Marks
- Basins are each 24"W x 20"L x 10"D with 30" W Basins available as an option
- 12" Backsplash
- (2) Pre-Rinse and Mixing Valve Faucet
- (1) Water Gun Kit (8 Tips)
- (1) Air Gun Kit (5 Tips)
- Deionized Faucet
- (2) Lever Handle Drain
- Optional Features include Overhead LED Lighting, 30"W Basins, Eye Wash, Mixing Valve Faucet, Pegboard Backframe.
- Limited 5 Year Warranty

SKU: LQDSA3012035

Categories: [decontamination-sinks](#)



PRODUCT DESCRIPTION

This 120" Wide Triple Basin Decontamination Sink provides ample room for cleaning. . With 300 Series Stainless Steel as a base, this Adjustable Height Decontamination Sink comes standard with a 12" Backsplash creating a universal workflow, 20"D x 24"W Basins, Pre-Rinse Faucet, 8-Tip Water Gun, 5-Tip Air Gun, Lever Handle Drain and more. And with a variety of optional accessories available including a Pegboard Backframe, 30" Basins, Overhead LED Lighting, these Decontamination Sinks can be customized to meet your SPD's needs.

ADDITIONAL INFORMATION

Dimensions	30 × 35 × 120 in
Department	sterile-processing
Product Line	Sinks
Product Category	decontamination-sinks
Product Subcategory	decontamination-sink-triple-basins
Material	STAINLESS-STEEL