



**DECONTAMINATION SINK, FIXED HEIGHT,  
96" WIDE, TRIPLE BASIN-  
(CAT.#LQDSF309637T)**

- 37" Fixed Height
- 96" Wide Triple Basin Sink with Etched Gallon Marks
- Basins are each 20"D x 24"W x 10"H with 30" W Basins available as an option
- 12" Backsplash
- (2) Pre-Rinse and Mixing Valve Faucet
- (1) Water Gun Kit (8 Tips)
- (1) Air Gun Kit (5 Tips)
- Deionized Faucet
- (2) Lever Handle Drain
- Optional Features include Overhead LED Lighting, 30"W Basins, Eye Wash, Mixing Valve Faucet, Pegboard Backframe.
- Limited 5 Year Warranty"

**SKU:** LQDSF309637T

**Categories:** [decontamination-sinks](#)

## PRODUCT DESCRIPTION

This 96" Wide Triple Basin Decontamination Sink packs Three Basins into a Double Basin Sink footprint, ideal for when space is at a premium. With 300 Series Stainless Steel as a base, this Fixed Height Decontamination Sink comes standard with a 12" Backsplash creating a universal workflow, 20"D x 24"W Basins, Pre-Rinse Faucet, 8-Tip Water Gun, 5-Tip Air Gun, Lever Handle Drain and more. And with a variety of optional accessories available including a Pegboard Backframe, 30" Basins, Overhead LED Lighting, these Decontamination Sinks can be customized to meet your SPD's needs.

## ADDITIONAL INFORMATION

**Dimensions**

30 x 96 x 37 in

**Department**

[sterile-processing](#)

**Product Line**

[Sinks](#)

**Product Category**

[decontamination-sinks](#)

