



## **DECONTAMINATION SINK, FIXED HEIGHT, DOUBLE BASIN-(CAT.#LQDSF309637D)**

- 37" Fixed Height
- Double Basin
- Double Basins with Etched Gallon Marks
- Basins are each 20"D x 24"W x 10"H with 30" W Basins available as an option
- 12" Backsplash
- (2) Pre-Rinse and Mixing Valve Faucet
- (1) Water Gun Kit (8 Tips)
- (1) Air Gun Kit (5 Tips)
- Deionized Faucet
- (2) Lever Handle Drain
- Optional Features include Overhead LED Lighting, 30"W Basins, Eye Wash, Mixing Valve Faucet, Pegboard Backframe.
- Limited 5 Year Warranty

---

**SKU:** LQDSF309637D

**Categories:** [decontamination-sinks](#)



### **PRODUCT DESCRIPTION**

The LogiQuip 300 Series Stainless Steel Fixed Height Decontamination Sinks comes standard with a 12" Backsplash creating a universal workflow, 20"D x 24"W Basins, Pre-Rinse Faucet, 8-Tip Water Gun, 5-Tip Air Gun, Lever Handle Drain and more. And with a variety of optional accessories available including a Pegboard Backframe, 30" Basins, Overhead LED Lighting, these Decontamination Sinks can be customized to meet your SPD's needs.

### **ADDITIONAL INFORMATION**

<b>Dimensions</b>	30 × 37 × 96 in
<b>Department</b>	<a href="#">sterile-processing</a>
<b>Product Line</b>	<a href="#">Sinks</a>
<b>Product Category</b>	<a href="#">decontamination-sinks</a>
<b>Product Subcategory</b>	<a href="#">decontamination-sink-double-basins</a>
<b>Material</b>	<a href="#">STAINLESS-STEEL</a>