



**DRI-STOR ENDOSCOPE CABINET, 16
LARGE DIAMETER SCOPE CAPACITY, BAR
CODE ACCESS-(CAT.#42590SCHG-16PTD-
TSBC)**

- Barcode Access with Endotrack
- HEPA Power Vent with Filter
- Channel Purges pushes air through each scope
- 304 Stainless Steel Construction
- 16 Durable, Easy-to-Clean Large Diameter Scope Holders
- LED Lighting
- Removable Drip Tray
- Limited 5-Year Warranty

SKU: 42590SCHG-16PTD-TSBC

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION



For highly efficient scope and probe drying, choose the Dri-Stor Drying Cabinet with HEPA Power Vent which circulates filtered air throughout the cabinet and through each scope for optimal storage conditions. Maximizes capacity for up to 16 large diameter scopes. Monitor and access equipment leveraging our Endotrack System for controlled access and scope tracking. Smooth-coated, stainless steel scope holders prevent damage upon hanging. Drip tray and holders can be removed easily for cleaning compliance.

ADDITIONAL INFORMATION

Dimensions	20 × 36 × 95 in
Department	cath-lab-ir-endoscopy
Product Line	Scope Cabinets
Product Category	dri-stor-scope-cabinets
Product Subcategory	dri-stor-16-scopes