



DRI-STOR ENDOSCOPE CABINET | 20 LARGE DIAMETER SCOPE CAPACITY- (CAT.#43890SCHG-20PSD)

- HEPA Power Vent and Channel Purge with filters
- 304 Stainless Steel Construction
- 20 Durable, Easy-to-Clean Large Diameter Scope Holders
- LED Lighting
- Removable Drip Tray
- Limited 5-Year Warranty

SKU: 43890SCHG-20PSD

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION

For best-in-class storage capacity and highly efficient drying, choose the Dri-Stor Ultra Endoscope Cabinet with HEPA Power Vent and Channel Purge. Maximizes capacity for up to 20 large diameter scopes. Monitor and access equipment through the tempered glass door. White acrylic scope holder for easy cleaning. Drip tray and acrylic scope holders can be removed easily for cleaning compliance.

ADDITIONAL INFORMATION

Dimensions	24 x 57 x 92 in
Department	cath-lab-ir-endoscopy
Product Line	Scope Cabinets
Product Category	dri-stor-scope-cabinets
Product Subcategory	dri-stor-20-scopes
Material	STAINLESS-STEEL
Construction	WELDED