



**DRI-STOR ENDOSCOPE CABINET | 20
LARGE DIAMETER SCOPE CAPACITY| RFID
ACCESS-(CAT.#43890SCHG-20PSD-TSRF)**

- RFID Access • HEPA Power Vent with Filter • 304 Stainless Steel Construction • 20 Durable, Easy-to-Clean Large Diameter Scope Holders • LED Lighting • Removable Drip Tray • Limited 5-Year Warranty

SKU: 43890SCHG-20PSD-TSRF

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION

For highly efficient drying, choose the Dri-Stor Endoscope Cabinet with HEPA Power Vent which circulates filtered air throughout the cabinet for optimal storage conditions. Maximizes capacity for up to 20 large diameter scopes. Monitor and access equipment through the tempered glass door. Smooth-coated, stainless steel scope holders prevent damage upon hanging. Drip tray and holders can be removed easily for cleaning compliance. Secure scopes with RFID access.

ADDITIONAL INFORMATION

Dimensions	24 × 57 × 92 in
Department	cath-lab-ir-endoscopy
Product Line	Scope Cabinets
Product Category	dri-stor-scope-cabinets
Product Subcategory	dri-stor-20-scopes
Material	STAINLESS-STEEL