



Product Image
Coming Soon

DRI-STOR ENDOSCOPE CABINET | 3 TEE PROBE CAPACITY-(CAT.#42380SCRG-3T)

- HEPA Power Vent with Filter • 304 Stainless Steel Construction • 3 Durable, Easy-to-Clean TEE Scope Holders • LED Lighting • Removable Drip Tray • Limited 5-Year Warranty

SKU: 42380SCRG-3T

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION

Up to 3 TEE Scopes dry and store inside this Dri-Stor Endoscope Cabinet with Overhead HEPA Power Vent. Designed to utilize vertical space of your suite, increasing endoscope capacity and drying efficiency. Monitor and access equipment through the tempered glass door. Smooth-coated, stainless steel scope holders prevent damage upon hanging. Drip tray and holders can be removed easily for cleaning compliance.

ADDITIONAL INFORMATION

Dimensions	18 × 27 × 87 in
Department	cath-lab-ir-endoscopy
Product Line	Scope Cabinets
Product Category	dri-stor-scope-cabinets
Product Subcategory	dri-stor-tee-probes
Material	STAINLESS-STEEL
Construction	WELDED