



DRI-STOR ENDOSCOPE CABINET | 6 SMALL DIAMETER SCOPE CAPACITY- (CAT.#42380SCRG-6PTD)

- HEPA Power Vent and Channel Purge with filters
- 304 Stainless Steel Construction
- 6 Durable, Easy-to-Clean Small Diameter Scope Holders
- LED Lighting
- Removable Drip Tray
- Limited 5-Year Warranty

SKU: 42380SCRG-6PTD

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION

Up to 6 small diameter scopes dry and store inside this Dri-Stor Endoscope Cabinet with Overhead HEPA Power Vent and Channel Purge. Designed to utilize vertical space of your suite, increasing endoscope capacity and drying efficiency. White acrylic scope holder for easy cleaning. Drip tray and acrylic scope holders can be removed easily for cleaning compliance.

ADDITIONAL INFORMATION

Dimensions	18 x 27 x 87 in
Department	cath-lab-ir-endoscopy
Product Line	Scope Cabinets
Product Category	dri-stor-scope-cabinets
Product Subcategory	dri-stor-6-scopes
Material	STAINLESS-STEEL
Construction	WELDED