



DRI-STOR FILTER KIT FOR HEPA-ONLY CABIENTS-(CAT.#40000HEPAKIT)

- Comes with 2 HEPA Pre-Filters, 2 HEPA Post-Filters, 2 O-Rings, 2 40000EF (External Filter)
- Replace every 6 onths

SKU: 40000HEPAKIT

Categories: [dri-stor-accessories](#)

PRODUCT DESCRIPTION

Keep your Dri-Stor Cabinet running like new with replacement filter kits. This Filter Kit comes with 2 HEPA Fan Pre-Filters, 2 HEPA Post-Filters, 2 O-Rings, and 2 External Filters.

ADDITIONAL INFORMATION

Department

[cath-lab-ir-endoscopy](#), [equipment-planners-architects](#), [new-construction-and-expansion](#), [sterile-processing](#)

Product Line

[Dri-Stor](#)

Product Category

[dri-stor-accessories](#)