



DRI-STOR SUPER ENDOSCOPE CABINET | 10 LARGE DIAMETER SCOPE CAPACITY- (CAT.#43690SCHG-10PS)

- HEPA Power Vent with Filter • 304 Stainless Steel Construction • 10 Durable, Easy-to-Clean Large Diameter Scope Holders • LED Lighting • Removable Drip Tray • Limited 5-Year Warranty


SKU: 43690SCHG-10PS

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION

Up to 10 large diameter scopes dry and store inside this Dri-Stor Endoscope Cabinet with HEPA Power Vent. Monitor and access equipment through the tempered glass door. Smooth-coated, stainless steel scope holders prevent damage upon hanging. Drip tray and holders can be removed easily for cleaning compliance.

ADDITIONAL INFORMATION



Dimensions	24 × 43 × 92 in
Department	cath-lab-ir-endoscopy
Product Line	Scope Cabinets
Product Category	dri-stor-scope-cabinets
Product Subcategory	dri-stor-10-scopes
Material	STAINLESS-STEEL
Construction	WELDED