



## **DRI-STOR ULTRA ENDOSCOPE CABINET | 20 LARGE DIAMETER SCOPE CAPACITY- (CAT.#43890SCHG-20PS)**

- HEPA Power Vent with Filter • 304 Stainless Steel Construction • 20 Durable, Easy-to-Clean Large Diameter Scope Holders • LED Lighting • Removable Drip Tray • Limited 5-Year Warranty

---

**SKU:** 43890SCHG-20PS

**Categories:** [dri-stor-scope-cabinets](#)

### **PRODUCT DESCRIPTION**

For best-in-class storage capacity and highly efficient drying, choose the Dri-Stor Ultra Endoscope Cabinet with HEPA Power Vent. Maximizes capacity for up to 20 large diameter scopes. Monitor and access equipment through the tempered glass door. Smooth-coated, stainless steel scope holders prevent damage upon hanging. Drip tray and holders can be removed easily for cleaning compliance.

### **ADDITIONAL INFORMATION**

<b>Dimensions</b>	24 x 57 x 92 in
<b>Department</b>	<a href="#">cath-lab-ir-endoscopy</a>
<b>Product Line</b>	<a href="#">Scope Cabinets</a>
<b>Product Category</b>	<a href="#">dri-stor-scope-cabinets</a>
<b>Product Subcategory</b>	<a href="#">dri-stor-20-scopes</a>
<b>Material</b>	<a href="#">STAINLESS-STEEL</a>
<b>Construction</b>	<a href="#">WELDED</a>

