



IV POLE, WEIGHTED BASE, 6-LEGS- (CAT.#IV46HW)

- 6-Leg Weighted Base • 4 Hooks, Stainless steel
- Adjusts from 52"-94" High • 3" Casters • White Base • Also available in Yellow, Green, Charcoal, Blue • Product Weight 22 lbs. • Limited 5-Year Warranty

SKU: IV46HW

Categories: [iv-poles](#)

PRODUCT DESCRIPTION

LogiQuip's Stainless steel Adjustable IV Poles have a weighted base that is 22" in diameter and 3" casters for stability and easy maneuverability allowing patients to move without interruption to intravenous therapies.

ADDITIONAL INFORMATION

| | |
|----------------------------|---|
| Dimensions | 22 x 5294 in |
| Department | emergency , patient-floor , surgical-services |
| Product Line | Complementary Products , IV Poles |
| Product Category | iv-poles |
| Product Subcategory | iv-poles-weighted-base |
| Material | STAINLESS-STEEL |
| Finish | STAINLESS-STEEL; POWDER-COATED |