



KICK BUCKET-(CAT.#KB)

With dolly. Removable 12 Quart basin. Perimeter bumper, Tip Resistant

SKU: KB

Categories: [kick-buckets](#)

PRODUCT DESCRIPTION

Stainless Steel, all-welded construction Kick Buckets are tip resistant and maneuver with ease on four (4) 2-inch Casters. These 12 quart capacity Kick Buckets with removable Buckets and Wrap-Around Bumper are an essential addition to any Operating Room, and available in two sizes.

ADDITIONAL INFORMATION

Dimensions	16 × 16 × 14 in
Department	emergency , sterile-processing
Product Line	Complementary Products
Product Category	kick-buckets
Material	STAINLESS-STEEL
Finish	STAINLESS-STEEL
Construction	ALL-WELDED

