



### **LOGICELL LOCKABLE DRAWER- (CAT.#2416LDR)**

- Locking Drawer
- Occupies 3 Cells
- Compatible with 16" Wide Column

---

**SKU:** 2416LDR

**Categories:** [logiccell-accessories](#), [parstor-accessories](#)

**Tags:** [LOGICELL](#)

## **PRODUCT DESCRIPTION**

Lockable Drawers are on full-extension slides and are used to store supplies in LogiCell Units . A top to enclose the drawer and a surface that engages with the lock are included.

## **ADDITIONAL INFORMATION**

### **Dimensions**

22 x 16 x 3.85 in

### **Department**

[cath-lab-ir-endoscopy](#), [sterile-processing](#), [supply-chain-clean-utility](#), [surgical-services](#)

### **Product Line**

[LogiCell](#)

### **Product Category**

[logiccell-accessories](#), [parstor-accessories](#)

### **Product Subcategory**

[logiccell-parstor-drawers](#)

### **Material**

[3/4 IN.; LAMINATE](#)

### **Finish**

[ALL-WELDED; STEEL; MELAMINE](#)

### **Color**

[FROSTY WHITE](#)

### **Tags**

[LOGICELL](#)