



MAYO STAND-(CAT.#M1621F)

- Stainless Steel Construction
- Singel Column Support
- 2" Casters
- Stainless Steel Removable Tray
- Friction Lock
- Limited 5-Year Warranty


SKU: M1621F

Categories: [mayo-stands](#)

PRODUCT DESCRIPTION

Keep surgical instruments and materials within reach during surgical procedures with LogiQuip Mayo Stands.

ADDITIONAL INFORMATION



Dimensions	16.25 × 21.25 × 58 in
Department	emergency , surgical-services
Product Line	Complementary Products
Product Category	mayo-stands
Product Subcategory	mayo-stands
Material	STEEL
Construction	FULLY WELDED