



MAYO STAND-(CAT.#MD1621F)


- Stainless Steel Construction
- Double Column Support
- Foot Operated
- 2" Casters
- Stainless Steel Removable Tray
- Friction Lock
- Limited 5-Year Warranty

SKU: MD1621F

Categories: [mayo-stands](#)

PRODUCT DESCRIPTION

Keep surgical instruments and materials within reach during surgical procedures with LogiQuip Mayo Stands.



ADDITIONAL INFORMATION

Dimensions	16.25 × 21.25 × 62 in
Department	emergency , surgical-services
Product Line	Complementary Products
Product Category	mayo-stands
Product Subcategory	mayo-stands
Material	STEEL
Construction	FULLY WELDED