



PARCART, 46 SLOTS, IV CART- (CAT.#500IV-V)

Kit Includes:

- (1) 5" Wire Basket with (2) Long Dividers
- (4) 8" Wire Baskets with (8) Long Dividers
- (1) Top Basket
- (4) 5" Casters-2 Swivel/2 Brake
- Holds Up To 100 Liter Bags
- 46 Slots
- Welded Solid Dolly Base with Open Frame
- Limited 5-Year Warranty

SKU: 500IV-V

Categories: [parcarts](#)

PRODUCT DESCRIPTION

An integrated Dolly Base with 5" Casters on an Open ParCarts creates stability and allows for the smooth transport of sterile supplies across your facility. This IV ParCart comes with Wire Baskets that are mounted on Full-Extension Slides that extend 100%, giving you full access to all of your supplies.

ADDITIONAL INFORMATION

Dimensions	27.5 × 25.5 × 58.5 in
Department	cath-lab-ir-endoscopy , sterile-processing , supply-chain-clean-utility , surgical-services
Product Line	ParStor
Product Category	parcarts
Product Subcategory	enclosed-carts
Material	STEEL
Finish	POWDER COATED
Construction	TUBULAR