



Product Image
Coming Soon

PARCART FRAME, 46 SLOTS, ENCLOSED WITH KEY LOCK-(CAT.#500HEN-KL)

- Enclosed ParCart with Key Lock • 46 Slots • Welded Solid Dolly Base • Empty • (4) 5" Casters-2 Swivel/2 Brake • Limited 5-Year Warranty

SKU: 500HEN-KL

Categories: [parcarts](#)

PRODUCT DESCRIPTION

An integrated Dolly Base with 5" Casters on an Enclosed ParCarts creates stability and allows for the smooth transport of sterile supplies across your facility. This secure ParCart can be customize with a variety of accessories mounted with Full-Extension Slides that extend 100%, giving you full access to all of your supplies.

ADDITIONAL INFORMATION

Dimensions	29.5 × 25.5 × 58.5 in
Department	cath-lab-ir-endoscopy , sterile-processing , supply-chain-clean-utility , surgical-services
Product Line	ParStor
Product Category	parcarts
Product Subcategory	enclosed-carts
Material	STEEL