



## **PARSTOR SHELVING KIT, 8-TIER, NARROW-(CAT.#NPS68CSM)**

Kit Includes: • Frame • (7) Waterfall Edge Full-Extension Shelves • (1) Waterfall Edge Top Shelf • (4) Casters, (2) Swivel, (2) Brake • (50) LH3, 3" Label Holders • Shelf Liner • 57 Slots • Add-On Unit Available, NPS68CSMA (Sold Separately) • Limited 5-Year Warranty

---

**SKU:** NPS68CSM

**Categories:** [parstor-carts](#)

## **PRODUCT DESCRIPTION**

This Narrow ParStor Starter Kit is ideal for use in Sterile Processing Departments (SPD) and Operating Rooms (OR). Our highest-density storage system, ParStor is designed to be used where space is at a premium. With a variety of storage accessories that extend 100%, you will have full access to surgical instruments, peel pouches, sutures, and any other needed supplies. Add-On units are available, sold separately.

## **ADDITIONAL INFORMATION**

<b>Dimensions</b>	24.5 × 18 × 74 in
<b>Department</b>	<a href="#">cath-lab-ir-endoscopy</a> , <a href="#">sterile-processing</a> , <a href="#">supply-chain-clean-utility</a> , <a href="#">surgical-services</a>
<b>Product Line</b>	<a href="#">ParStor</a>
<b>Product Category</b>	<a href="#">parstor-carts</a>
<b>Product Subcategory</b>	<a href="#">open-storage</a>