



RING STAND-(CAT.#DBS)

- (2) 7-Quart Stainless Steel Basins
- (4) 2" Casters
- H-Brace Support
- Limited 5-Year Warranty

SKU: DBS

Categories: [solution-stand](#)

PRODUCT DESCRIPTION

Stainless Steel, fully-welded Double Basin Solution Stand holds two 7 quart basin. Used in hospitals and surgery centers, these Ring Stands come with an H-Brace for added stability.

ADDITIONAL INFORMATION

Dimensions	14 × 28 × 34 in
Department	emergency, surgical-services
Product Line	Complementary Products
Product Category	solution-stand
Product Subcategory	solution-stand
Material	STEEL
Construction	FULLY WELDED