



SCOPE CABINET-(CAT.#43490SCHG-10PT)

- HEPA Power Vent with Filter • SGNA and AORN Compliant • 304 Stainless Steel Construction • 10 Durable, Easy-to-Clean Large Diameter Scope Holders • LED Lighting • Removable Drip Tray • Maintenance Plans Available • Limited 5-Year Warranty


SKU: 43490SCHG-10PT

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION

Dri-Stor Top Mount Scope Cabinets are equipped with a HEPA Power Vent circulates HEPA Filtered Air throughout the Cabinet for infection-control conscious drying and storage The resulting pressurized environment creates a barrier, preventing non-filtered air from entering

ADDITIONAL INFORMATION



Dimensions	24 × 30 × 92 in
Department	cath-lab-ir-endoscopy
Product Line	Scope Cabinets
Product Category	dri-stor-scope-cabinets
Product Subcategory	dri-stor-10-scopes
Material	STAINLESS-STEEL
Construction	WELDED