



SCOPE CABINET- **(CAT.#43690SCHG-20PTD)**

- HEPA Power Vent with Channel Purge •
- ANSI/AAMI ST91: 2021, SGNA, and AORN
- Compliant • 304 Stainless Steel Construction •
- 20 Durable, Easy-to-Clean Large Diameter Scope
- Holders • LED Lighting • Removable Drip Tray •
- Maintenance Plans Available • Limited 5-Year
- Warranty

SKU: 43690SCHG-20PTD

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION

HEPA Power Vent with Channel Purge Technology comply with ANSI/AAMI ST91:2021 Guidelines, emphasizing the importance of drying the internal channels of flexible and semi-rigid endoscopes by forcing HEPA Filtered Air through the tubing channels assigned to individual Scopes for 30-minute cycles Airflow is discontinued for 20 minutes before a new cycle begins The Filtered Air quality is to 0.3 micron at 99.97% efficiency.

ADDITIONAL INFORMATION

Dimensions	24 × 44 × 92 in
Department	cath-lab-ir-endoscopy
Product Line	Scope Cabinets
Product Category	dri-stor-scope-cabinets
Product Subcategory	dri-stor-20-scopes
Material	STAINLESS-STEEL
Construction	WELDED