



SCOPE CABINET, STAINLESS STEEL, TOP MOUNT, 10 LARGE DIAMETER SCOPE CAPACITY-(CAT.#42490SCHG-10PT)

- HEPA Power Vent with Filter
- SGNA and AORN Compliant
- 304 Stainless Steel Construction
- 10 Durable, Easy-to-Clean Large Diameter Scope Holders
- LED Lighting • Removable Drip Tray
- Maintenance Plans Available
- Limited 5-Year Warranty

SKU: 42490SCHG-10PT

Categories: [dri-stor-scope-cabinets](#)

PRODUCT DESCRIPTION

Dri-Stor Top Mount Scope Cabinets are equipped with a HEPA Power Vent circulates HEPA Filtered Air throughout the Cabinet for infection-control conscious drying and storage. The resulting pressurized environment creates a barrier, preventing non-filtered air from entering

ADDITIONAL INFORMATION

Department	cath-lab-ir-endoscopy
Product Line	Scope Cabinets
Product Category	dri-stor-scope-cabinets
Product Subcategory	dri-stor-10-scopes
Material	STAINLESS-STEEL