



#### **TABLE-(CAT.#LQBT1630)**

- Fully Welded Construction
- 16 Gauge, #304 Stainless Steel
- Reinforced Tabletops
- 1" Square Posts
- 3" Casters
- Lower Shelf for Additional Storage and Support
- Limited 5-Year Warranty

---

**SKU:** LQBT1630

**Categories:** [back-instrument-stainless-tables](#)

## **PRODUCT DESCRIPTION**

LogiQuip Worktables (also known as Instrument Tables and Back Tables) come in several sizes. This 16" x 30" x 34" Stainless Steel Back Table has fully welded construction and is ideal for sterile processing, surgical services, pharmacy, and cath lab usage.



## **ADDITIONAL INFORMATION**

<b>Dimensions</b>	16 × 30 × 34 in
<b>Department</b>	<a href="#">cath-lab-ir-endoscopy</a> , <a href="#">emergency</a> , <a href="#">surgical-services</a>
<b>Product Line</b>	<a href="#">Complementary Products</a>
<b>Product Category</b>	<a href="#">back-instrument-stainless-tables</a>
<b>Product Subcategory</b>	<a href="#">back-instrument-stainless-tables</a>
<b>Material</b>	<a href="#">STAINLESS-STEEL</a>
<b>Finish</b>	<a href="#">STAINLESS-STEEL</a>
<b>Color</b>	<a href="#">STAINLESS-STEEL</a>