



TABLE-(CAT.#LQBT2436)

- Fully Welded Construction
- 16 Gauge, #304 Stainless Steel
- Reinforced Tabletops
- 1" Square Posts
- 3" Casters
- Lower Shelf for Additional Storage and Support
- LQBT2448 image shown
- Limited 5-Year Warranty

SKU: LQBT2436

Categories: [back-instrument-stainless-tables](#)

PRODUCT DESCRIPTION

LogiQuip Worktables (also known as Instrument Tables and Back Tables) come in several sizes. This 24" x 36" x 34" Stainless Steel Back Table has fully welded construction and is ideal for sterile processing, surgical services, pharmacy, and cath lab usage.

ADDITIONAL INFORMATION

Dimensions	24 × 36 × 34 in
Department	cath-lab-ir-endoscopy , emergency , surgical-services
Product Line	Complementary Products
Product Category	back-instrument-stainless-tables
Product Subcategory	back-instrument-stainless-tables
Material	STAINLESS-STEEL
Finish	STAINLESS-STEEL
Color	STAINLESS-STEEL