



TREATMENT CART, 36.5" TALL NARROW
TREATMENT CART, KEY LOCK, FOUR
SOFT-CLOSE DRAWERS-
(CAT.#QTN37W24K3669W)

- 4 Soft-Close Drawers, (1) 3", (2) 6", and (1) 9" -
- Includes (1) 3" and (1) 6" Drawer Divider Set •
- Pull-Out Shelf For Added Workspace • eLock for
- Secure Keyless Entry • 5" Hospital Grade Casters
- (2 Swivel, 1 Brake, 1 Directional Lock) • QTN40
- Image Shown • Several Optional Accessories
- Available including Glove Dispenser, Waste
- Basket, and Sharps Container Set • Limited 5
- Year Warranty

SKU: QTN37W24K3669W

Categories: [qt-carts](#)

PRODUCT DESCRIPTION

Designed for general use, this Double Wide Quiet Treatment Procedure Cart comes standard with Key Lock for secure storage, eight Soft-Closing Drawers and Hospital Grade Casters to help provide a quieter environment for patient care.

ADDITIONAL INFORMATION

Dimensions	21 × 20 × 36.5 in
Department	cath-lab-ir-endoscopy , emergency , equipment-planners-architects , linen-environmental-services , new-construction-and-expansion , patient-floor , sterile-processing , supply-chain-clean-utility , surgical-services
Product Line	QT Carts
Product Category	QT Carts
Product Subcategory	enclosed-carts
Material	ALUMINUM