



TREATMENT CART, 41" TALL, KEY LOCK, FIVE SOFT-CLOSE DRAWERS, PULL-OUT SHELF-(CAT.#QT41W27K33669W)

- 5 Soft-Close Drawers, (2) 3", (2) 6", and (1) 9" - Includes (1) 3" and (1) 6" Drawer Divider Set • Pull-Out Shelf For Added Workspace • eLock for Secure Keyless Entry • 5" Hospital Grade Casters (2 Swivel, 1 Brake, 1 Directional Lock) • Several Optional Accessories Available including Glove Dispenser, Waste Basket, and Sharps Container Set • Limited 5 Year Warranty

SKU: QT41W27K33669W

Categories: [qt-carts](#)

PRODUCT DESCRIPTION

Designed for general use, this Quiet Treatment Procedure Cart come standard with Key Lock for secure storage, five Soft-Closing Drawers and Hospital Grade Casters to help provide a quieter environment for patient care.

ADDITIONAL INFORMATION

Dimensions	20 × 27 × 40.5 in
Department	cath-lab-ir-endoscopy , emergency , equipment-planners-architects , linen-environmental-services , new-construction-and-expansion , patient-floor , sterile-processing , supply-chain-clean-utility , surgical-services
Product Line	QT Carts
Product Category	QT Carts
Product Subcategory	enclosed-carts
Material	ALUMINUM