



**TREATMENT CART, 43" TALL NARROW,
KEY LOCK, SIX SOFT-CLOSE DRAWERS,
PULL-OUT SHELF-
(CAT.#QTN43W30K333669W)**

- 4 Soft-Close Drawers, (3) 3", (2) 6", and (1) 9" - Includes (1) 3" and (1) 6" Drawer Divider Set • Pull-Out Shelf For Added Workspace • eLock for Secure Keyless Entry • 5" Hospital Grade Casters (2 Swivel, 1 Brake, 1 Directional Lock) • QTN40 Image Shown • Several Optional Accessories Available including Glove Dispenser, Waste Basket, and Sharps Container Set • Limited 5 Year Warranty

SKU: QTN43W30K333669W

Categories: [qt-carts](#)

PRODUCT DESCRIPTION

Designed for general use, this Double Wide Quiet Treatment Procedure Cart comes standard with Key Lock for secure storage, eight Soft-Closing Drawers and Hospital Grade Casters to help provide a quieter environment for patient care.

ADDITIONAL INFORMATION

Dimensions	21 × 20 × 42.5 in
Department	cath-lab-ir-endoscopy , emergency , equipment-planners-architects , linen-environmental-services , new-construction-and-expansion , patient-floor , sterile-processing , supply-chain-clean-utility , surgical-services
Product Line	QT Carts
Product Category	QT Carts
Product Subcategory	enclosed-carts
Material	ALUMINUM