



**TREATMENT CART, 44" TALL PEDIATRIC
TREATMENT CART, INDIVIDUAL
BREAKAWAY LOCKS, NINE SOFT-CLOSE
DRAWERS, PULL-OUT SHELF-
(CAT.#QT44W30BPEDS2)**

- 9 Soft-Close Drawers, (8) 3" and (1) 6" - Includes (1) 3" and (1) 6" Drawer Divider Set • Pull-Out Shelf For Added Workspace • Individual Breakaway Locks (one for each drawer) with Key Lock • 5" Hospital Grade Casters (2 Swivel, 1 Brake, 1 Directional Lock) • Several Optional Accessories Available including Glove Dispenser, Waste Basket, and Sharps Container Set • Limited 5 Year Warranty

SKU: QT44W30BPEDS2

Categories: [qt-carts](#)

PRODUCT DESCRIPTION

Designed for use in Pediatrics, this Quiet Treatment Pediatric Treatment Cart comes standard with individual breakaway locks (one for each drawer) for quick access to secure items, nine Soft-Closing Drawers and Hospital Grade Casters to help provide a quieter environment for patient care.

ADDITIONAL INFORMATION

Dimensions	20 × 27 × 43.5 in
Department	cath-lab-ir-endoscopy , emergency , equipment-planners-architects , linen-environmental-services , new-construction-and-expansion , patient-floor , sterile-processing , supply-chain-clean-utility , surgical-services
Product Line	QT Carts
Product Category	QT Carts
Product Subcategory	enclosed-carts
Material	ALUMINUM