



## **WARMING CABINET-(CAT.#W2624TD)**

W-Series Data Logging Easy Log Integrated Data Logger with removable USB-standard! Automatically records data every 30 minutes. Internal storage up to 2 years. Digitally Controlled Heating Chambers. Locking temperature set point. Heat control accuracy to plus or minus 3 degrees F. Rapid warm time (2-6 hours depending on load). Uniform heating. Heating range 90-160 degrees F. Fully Insulated. U.L./C.U.L. approved. Audible and visual over/under temp alarm. Compartment size: 13.24 cubic ft. Electrical: 120V, 12 Amp., 50/60HZ. FDA registered and OSHPD pre-approved.

---

**SKU:** W2624TD

**Categories:** [warmers](#)



## **PRODUCT DESCRIPTION**

Double Door Compartment, Solid stainless-steel Doors, 3 adjustable shelves. (Pictured with glass doors.)

## **ADDITIONAL INFORMATION**

|                            |                                            |
|----------------------------|--------------------------------------------|
| <b>Dimensions</b>          | 26.5 × 24 × 74.5 in                        |
| <b>Department</b>          | <a href="#">supply-chain-clean-utility</a> |
| <b>Product Line</b>        | <a href="#">Complementary Products</a>     |
| <b>Product Category</b>    | <a href="#">warmers</a>                    |
| <b>Product Subcategory</b> | <a href="#">warmers</a>                    |
| <b>Material</b>            | <a href="#">STAINLESS-STEEL</a>            |
| <b>Finish</b>              | <a href="#">STAINLESS-STEEL</a>            |