



WARMING CABINET-(CAT.#W2630MDG)

W-Series Data Logging Easy Log Integrated Data Logger with removable USB-standard!
Automatically records data every 30 minutes.
Internal storage up to 2 years. Digitally Controlled Heating Chambers. Locking temperature set point. Heat control accuracy to plus or minus 3 degrees F. Rapid warm time (2-6 hours depending on load). Uniform heating. Heating range 90-160 degrees F. Fully Insulated. U.L./C.U.L. approved. Audible and visual over/under temp alarm. Compartment size: 9.72 cubic ft. Electrical: 120V, 12 Amp., 50/60HZ. FDA registered and OSHPD pre-approved. Weight 484 lbs.

SKU: W2630MDG

Categories: [warmers](#)



PRODUCT DESCRIPTION

Double Door Compartment, Glass Doors, 3 adjustable shelves. (Pictured with optional dolly.)

ADDITIONAL INFORMATION

Dimensions	20.5 × 30 × 64.75 in
Department	supply-chain-clean-utility
Product Line	Complementary Products
Product Category	warmers
Product Subcategory	warmers
Material	STAINLESS-STEEL
Finish	STAINLESS-STEEL