



WARMING CABINET-(CAT.#W2630TT)

W-Series Data Logging Easy Log Integrated Data Logger with removable USB-standard! Automatically records data every 30 minutes. Internal storage up to 2 years. Digitally Controlled Heating Chambers. Locking temperature set point. Heat control accuracy to plus or minus 1 degrees F. Rapid warm time (2-6 hours depending on load). Uniform heating. Heating range: top 90-110, middle 90-135, bottom 90-160 degrees F. Fully Insulated. U.L./C.U.L. approved. Audible and visual over/under temp alarm. Compartment size: 14.01 cubic ft. Electrical: 120V, 13.5 Amp., 50/60HZ. FDA registered and OSHPD pre-approved.

SKU: W2630TT

Categories: [warmers](#)



PRODUCT DESCRIPTION

Full size, triple full-view glass doors. 1 adjustable shelf.

ADDITIONAL INFORMATION

Dimensions	26.5 × 30 × 78.75 in
Department	supply-chain-clean-utility
Product Line	Complementary Products
Product Category	warmers
Product Subcategory	warmers
Material	STAINLESS-STEEL
Finish	STAINLESS-STEEL