



WARMING CABINET-(CAT.#W2630XSS)

W-Series Data Logging Easy Log Integrated Data Logger with removable USB-standard!
Automatically records data every 30 minutes. Internal storage up to 2 years. Digitally Controlled Heating Chambers. Locking temperature set point. Heat control accuracy to plus or minus 1 degree F. Rapid warm time (2-6 hours depending on load). Uniform heating. Heating range 90-160 degrees F. Fully Insulated. U.L./C.U.L. approved. Audible and visual over/under temp alarm. Compartment size: 5.27 cubic ft. Electrical: 120V, 7.5 Amp., 50/60HZ. FDA registered and OSHPD pre-approved. Weight 202 lbs.

SKU: W2630XSS

Categories: [warmers](#)

PRODUCT DESCRIPTION



Small, single solid stainless-steel door, 1 adjustable shelf. (Pictured with glass doors and additional shelf)

ADDITIONAL INFORMATION

Dimensions	26.5 × 30 × 24.5 in
Department	supply-chain-clean-utility
Product Line	Complementary Products
Product Category	warmers
Product Subcategory	warmers
Material	STAINLESS-STEEL
Finish	STAINLESS-STEEL