

## **WIRE SHELF, 24" X 30"-(CAT.#2430CH)**



- Double Chrome Plated, Carbon Steel Core
- Resists discoloration and peeling
- Compatible with LQS 1" Round Posts
- Max Weight Capacity 800 lbs.
- NSF Certified
- Limited 5-Year Warranty


---

**SKU:** 2430CH

**Categories:** [wire-shelving-components](#)

### **PRODUCT DESCRIPTION**

LogiQuip's Standard (LQS) line of Flat Wire Shelving begins with a carbon steel core that is double chrome plated for strength and durability that lasts. Leverage our 1" Posts, LQS Flat Wire, Solid Stainless Steel, and Waterfall Edge Shelving as well as our collared Basket Shelving to create a custom shelving unit that can be used virtually anywhere within your facility.



### **ADDITIONAL INFORMATION**

<b>Dimensions</b>	24 x 30 in
<b>Department</b>	<a href="#">cath-lab-ir-endoscopy</a> , <a href="#">sterile-processing</a> , <a href="#">supply-chain-clean-utility</a> , <a href="#">surgical-services</a>
<b>Product Line</b>	<a href="#">Kanban Shelving</a> , <a href="#">Wire Shelving Carts</a>
<b>Product Category</b>	<a href="#">wire-shelving-components</a>
<b>Product Subcategory</b>	<a href="#">wire-shelving-shelves</a>
<b>Material</b>	<a href="#">STEEL</a>
<b>Finish</b>	<a href="#">CHROME</a>