



WORKSTATION, 72" W X 36" D ELECTRIC HEIGHT-ADJUSTABLE-(CAT.#HTE7236)

- Patented Electronic Control Height Adjustment Workstation • Height Range 34" - 52.5" • (4) 5" Total-Lock Casters • Powder Coated Steel Frame
- Stainless Steel Work Surface • Maximum Weight is 400 Pounds • Limited 5-Year Warranty

SKU: HTE7236

Categories: [workstation-height-adjustable](#)

PRODUCT DESCRIPTION

Easily raise or lower your work surface with LogiQuip's fully electric, adjustable Workstations. With the press of a button, the work surface can be adjusted from 34"-52.5", providing an ergonomic Workstation that allows the user to sit or stand while working. Ideal for Sterile Processing, these prep and pack Workstations feature a powder-coated steel frame, Stainless Steel work surface, and two slotted uprights and top frame that accommodate several optional accessories making this Workstation fully customizable.



ADDITIONAL INFORMATION

Dimensions	36 × 72 × 52.5 in
Department	equipment-planners-architects , new-construction-and-expansion , sterile-processing , supply-chain-clean-utility , surgical-services
Product Line	Workstations
Product Category	workstation-height-adjustable
Product Subcategory	workstations
Finish	PAINTED-METAL