



WRAP RACK-(CAT.#WC60CH)

Kit Includes:

- (4) Adjustable Bars with Rubber O Rings
- (4) 5" Casters - (2) Swivel and (2) Swivel with Brake
- (4) Donut Bumpers
- Three Sided Frame
- Bottom Shelf and Shelf Cover
- Double Chrome Plated, Carbon Steel Core Wire Shelves, Basket Shelves, and Posts
- Resists discoloration and peeling.
- Limited 5-Year Warranty

SKU: WC60CH

Categories: [wire-shelving-carts](#)

PRODUCT DESCRIPTION

LogiQuip's 24" x 60" Sterile Wrap Carts with (4) adjustable bars and rubber rings to keep wrap in place provide a convenient way to store and access sterile wrap.

ADDITIONAL INFORMATION

Dimensions	24 × 60 × 63 in
Department	sterile-processing
Product Line	Wire Shelving
Product Category	wire-shelving-carts
Product Subcategory	wire-racks-wrap-racks
Material	STEEL
Finish	CHROME