

## The Industry Demands Drying Cabinets!

The cost of an infection from an improperly cleaned/stored endoscope can be significant. Recent outcomes have driven professional organizations such as the SGNA and AORN to rewrite guidelines for endoscope reprocessing. Recognizing that drying is an element in reprocessing, they state that endoscopes must be stored in an area that is clean, well-ventilated and dust-free to keep them dry and free of microbial contamination. Furthermore, endoscopes should be stored in a closed cabinet with positive pressure that allows circulation of filtered air around the endoscopes when a cabinet with channel drying is not available.

## Supportive Design

LogiQuip® is proud to offer a full line of compliant, positive-pressure cabinets with HEPA filtered circulated air in economical melamine construction. Accommodating large and small diameter scopes and TEE probes, LogiQuip® cabinets hold 10% more scopes than most competitive models. An expansive offering includes clear doors for visual access, a choice of cabinet widths, lock types, and various types of device holders. Our Dri-Stor cabinets have a broad range of options including channel air purge. Whichever combination you select, we can help you protect and manage your scopes.

## Scope Cabinet Product Offering

MODEL #	DIMENSIONS	SCOPE	CAPACITY	DOOR
	W x H (inches)	diameter	scopes	
SCL-L27S-CL	27 x 92	large	6, single row	left hinge glass door
SCL-L27S-CR	27 x 92	large	6, single row	right hinge glass door
SCL-L27D-CL	27 x 92	large	12 (6 per row)	left hinge glass door
SCL-L27D-CR	27 x 92	large	12 (6 per row)	right hinge glass door
SCL-L36S-CD	36 x 92	large	9, single row	double glass doors
SCL-L36D-CD	36 x 92	large	18 (9 per row)	double glass doors
SCL-S27S-CL	27 x 60	small	5, single row	left hinge glass door
SCL-S27S-CR	27 x 60	small	5, single row	right hinge glass door
SCL-S27D-CL	27 x 60	small	10 (5 per row)	left hinge glass door
SCL-S27D-CR	27 x 60	small	10 (5 per row)	right hinge glass door
SCL-S36S-CD	36 x 60	small	7, single row	double glass doors
SCL-S36D-CD	36 x 60	small	14 (7 per row)	double glass doors
SCL-T27S-CL	27 x 84	tee probe	3	left hinge glass door
SCL-T27S-CR	27 x 84	tee probe	3	right hinge glass door
SCL-T36S-CD	36 x 84	tee probe	5	double glass doors



# Scope Drying Cabinets

## Scope Cabinet Product Offering Cont.



Cabinet for small diameter scopes

### Compliant Design

- Constant fan activity for maximum circulation and drying
- Premium HEPA filter removes 99.97% of particulates 0.3 micron
- Lockable doors
- Glass doors for faster visual access
- Optional sloped top to support infection control
- Removable stainless steel drip tray allows access to the floor for optimal cleaning
- Especially designed holder cradles and protects the scopes
- Padded rear wall protects the scopes during hanging and removal
- Cord retainer prevents scopes from touching each other

### Innovative Design

- Dri-Stor offers more ventilation options than any other cabinet made in North America
- Dri-Stor cabinets feature both positive pressure and individual channel purge drying
- Dri-Stor is easily upgradable, so you can start with basic cabinet and move up to channel purge later
- LED Lighting
- Two cabinets with channel purge can draw from a single component cabinet, providing significant space/cost savings

***Please contact your LogiQuip® rep for more information about our Dri-Stor product lines.***

LogiQuip® has been in the healthcare storage and distribution systems business since 1992. Our value-added services support the hospital clinician, logistics officer, purchaser and planner.

Because each hospital's inventory, supply, and storage requirements are diverse, our mission is to solve your specific storage and distribution issues. With "space" being one of the hospital's most expensive capital assets, our philosophy is to maximize storage efficiency while supporting infection control and functional ergonomics. Whether a single cart or a complete supply storage system, our experienced healthcare storage solutions adviser can offer a solution to your challenges.



Easy access makes filter replacement a breeze



www.logiquip.com info@logiquip.net

Phone: (800) 665-3760

Fax: (888) 965-3984